

### Care-recommendation

material: Skinfoam - coating



**Advantages:**

- > Polyurethan-coating which is kind to the skin
- > very tough surface
- > special versions are possible

**X-ray possibility:**

- > from 70 –140 KV
- > standard colours yellow, blue and green

**Colours:**

- > x-ray standard colour is yellow
- > care standard colour is skin colour
- > silver-grey is an antistatic coating
- > other colours please infer from our brochure

**Care:**

- > wipe damply

**Disinfection:** undiluted Antifect Liquid (Schülke & Mayr)

**ATTENTION:** The material is dampness-reject but not absolutely waterproof. Please avoid in matter of this fact wetness to prolong the life of your Polyform-product. We cannot guarantee a absolutely cleaning of some substances (e.g. blood, secretion etc.).