



IDEEN & TECHNIK  
IN SCHAUMSTOFF

### Care-recommendation

**material:** Skinfoam - coating  
(out of Polyurethan)



**advantages:**

- > Polyurethan-coating which is kind to the skin
- > tough surface
- > custom-made versions are possible

**x-ray possibility:**

- > from 70 –140 KV
- > standard colours yellow, blue and green

**colours:**

- > x-ray standard colour is yellow
- > care standard colour is skin colour
- > silver-grey is an antistatic coating
- > other colours please infer from our brochure

**care:**

- > wipe damply

**disinfection:** undiluted Antifect Liquid (Schülke & Mayr)

**ATTENTION:** The material is dampness-reject but not absolutely waterproof. Please avoid in matter of this fact wetness to prolong the life of your Polyform-product. We cannot guarantee a absolutely cleaning of some substances (e.g. blood, secretion etc.).