



IDEAS & TECHNIQUES IN FOAM

POLYFORM Kunststofftechnik
GmbH & Co. Betriebs KG



Positioning aids
Healthcare applications

Contents

04 - 05

LEG SPLINTS

06 - 07

SEATING CUSHIONS & HEEL SUPPORTS

08 - 11

POSITIONING ROLLS

12 - 13

HEAD SUPPORTS

14

POSITIONING AIDS &
MATTRESSES

15

ARM SUPPORTS

16

TRIANGLES

17

RECTANGLES

18

WEDGES

19

CIRCULAR FOAM PADS

20 - 22

POSITIONING-AIDS-SET

23

SINGLE-USE PRODUCTS

24

TRAINING ARTICLES

25

COVERING FABRICS & COATINGS
CUSTOM-MADE PRODUCTS

26

SKINFOAM

POSITIONING AIDS

Leg splints



Leg splint 30° - wide form

Dimensions	90 x 26 x 20 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115016	4001363000283
	Skinfoam blue	115018	4001363000290
	Skinfoam violet (ET176)	128593	4001363003208
Material Colour	foam + imitation leather		
	jade	125024	4001363002294
	black	125023	4001363002300
	ocean	127772	4001363002317
	white (ET176)	127773	4001363002324




Leg splint 45° - wide form

Dimensions	72 x 26 x 20 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115020	4001363000306
	Skinfoam blue	117975	4001363002614
	Skinfoam violet (ET176)	128594	4001363003215
Material Colour	foam + imitation leather		
	jade	125025	4001363002355
	black	125026	4001363002362
	ocean	127776	4001363002348
	white (ET176)	127775	4001363002331




Leg splint flat - wide form

Dimensions	74 x 30 x 14 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	112827	4001363002621
	Skinfoam blue	115045	4001363000337
	Skinfoam violet (ET176)	128596	4001363003222
Material Colour	foam + imitation leather		
	jade	125027	4001363002379
	black	125028	4001363002409
	ocean	107218	4001363002386
	white (ET176)	127777	4001363002393



Leg splint flat - narrow form

Dimensions	67 x 22 x 10,5 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115046	4001363000344
	Skinfoam blue	115047	4001363000351
	Skinfoam violet (ET176)	128597	4001363003239
Material Colour	foam + imitation leather		
	jade	125029	4001363002430
	black	125030	4001363002447
	ocean	127780	4001363002423
	white (ET176)	127779	4001363002416



Leg splint "pius"
for adults

Dimensions	67 x 27 x 13 cm	Item No.:	UDI
Material Colour	Skinfoam blue (ET176)	103797	4001363000979



Leg splint "pius"
for adults, longer version

Dimensions	80 x 29 x 13 cm	Item No.:	UDI
Material Colour	Skinfoam blue (ET176)	115268	4001363000412



Leg splint "pius"
for adults, shorter version

Dimensions	30 x 26,5 x 13 cm	Item No.:	UDI
Material Colour	Skinfoam blue (ET176)	115265	4001363000405



Leg splint "pius"
for children

Dimensions	53 x 22 x 11 cm	Item No.:	UDI
Material Colour	Skinfoam blue (ET176)	103798	4001363000986



Heel support

Dimensions	26,5 x 16 x 9 cm	Item No.:	UDI
Material Colour	Skinfoam blue (ET176)	128591	4001363003239



Trade purchase prices net per piece plus VAT

POSITIONING AIDS

Seating cushions



Seating cushion round

Dimensions	Ø 45 x 7 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115834 4001363000634 103796 4001363000993 128598 4001363003246	



Seating cushion oval

Dimensions	Ø 47 x 38 x 7 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115835 4001363000641 115836 4001363000658 128599 4001363003253	



Seating cushion U-shaped

Dimensions	40 x 38 x 7 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115838 4001363000665 115839 4001363000672 128600 4001363003260	



Expanding wedge for adults

Dimensions	38/15 x 30 x 20 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115048 4001363000368 115049 4001363000375 128601 4001363003277	



Expanding wedge for children


Dimensions	36/15 x 25 x 12 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115050 4001363000382 115052 4001363000399 128602 4001363003284	



Trade purchase prices net per piece plus VAT

Wedge for 30°- positioning

Dimensions	60 x 30 / 12 x 12 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115840 115841 128603	4001363000689 4001363000696 4001363003291



Wedge for 30°- positioning
with cotton cretonne cover

Dimensions	60 x 30 / 12 x 12 cm	Item No.:	UDI
Material	ET176	129846	4001363003901




Heel support with hole | pair

Dimensions	∅10 / 3 x 4 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115843 103813 128604	4001363000702 4001363001006 4001363003307



Heel support with hole | pair

Dimensions	∅15 / 5 x 4 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115844 115167 128605	4001363000719 4001363001013 4001363003314



Heel support with hole | pair

Dimensions	∅20 / 8 x 4 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115846 115168 128606	4001363000726 4001363001020 4001363003321



Trade purchase prices net per piece plus VAT

POSITIONING AIDS

Positioning rolls



Positioning roll

Dimensions	Ø10 x 25 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103891	4001363000054
Material Colour	Skinfoam yellow	103892	4001363000078
	Skinfoam blue	103893	4001363000061
	Skinfoam violet (ET176)	113016	4001363003185
Material Colour	foam + imitation leather		
	jade	125035	4001363002157
	black	125034	4001363002164
	ocean	127517	4001363002560
	white (ET176)	125033	4001363002140
Dimensions	Ø15 x 40 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	107471	4001363000085
	Skinfoam blue	115269	4001363000429
	Skinfoam violet (ET176)	128607	4001363003338
Material Colour	foam + imitation leather		
	jade	125037	4001363002195
	black	125038	4001363002201
	ocean	127754	4001363002171
	white (ET176)	125036	4001363002539



Trade purchase prices net per piece plus VAT

Positioning roll

Dimensions	Ø 15 x 60 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115372	4001363000436
	Skinfoam blue	115373	4001363000443
	Skinfoam violet (ET176)	128608	4001363003345
Material Colour	foam + imitation leather	jade	110702 4001363002584
		black	125040 4001363002607
		ocean	127889 4001363002591
		white	125039 4001363002577
		(ET176)	
Dimensions	Ø 22 x 40 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	103897	4001363000092
	Skinfoam blue	115374	4001363000450
	Skinfoam violet (ET176)	128609	4001363003352
Material Colour	foam + imitation leather	jade	125042 4001363002232
		black	125043 4001363002249
		ocean	127759 4001363002218
		white	125041 4001363002225
		(ET176)	
Dimensions	Ø 22 x 60 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	103898	4001363000108
	Skinfoam blue	115170	4001363002836
	Skinfoam violet (ET176)	128610	4001363003369
Material Colour	foam + imitation leather	jade	125046 4001363002270
		black	125047 4001363002287
		ocean	125044 4001363002256
		white	125045 4001363002263
		(ET176)	



Trade purchase prices net per piece plus VAT

POSITIONING AIDS

Positioning rolls



1/2 Positioning roll


Dimensions	Ø 15/7,5 x 40 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103900	4001363000115
Material Colour	Soft-PE violet	122092	4001363001549
Material Colour	Skinfoam yellow	110954	4001363000122
	Skinfoam blue	115375	4001363000467
	Skinfoam violet (ET176)	128611	4001363003376
Material Colour	foam + imitation leather		
	jade	127739	4001363002003
	black	125048	4001363002210
	ocean	127738	4001363001990
	white (ET176)	127736	4001363000139
Dimensions	Ø 22/11 x 40 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	107236	4001363000153
	Skinfoam blue	107235	4001363000146
	Skinfoam violet (ET176)	128612	4001363003383
Material Colour	foam + imitation leather		
	jade	127344	4001363002553
	black	127341	4001363002546
	ocean	127339	4001363000160
	white (ET176)	127340	4001363002539



Trade purchase prices net per piece plus VAT

3/4 Positioning roll

Dimensions	Ø 15/11 x 40 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115378	4001363000474
	Skinfoam blue	111944	4001363001037
	Skinfoam violet (ET176)	128613	4001363003390
Material Colour	foam + imitation leather		
	jade	125051	4001363002034
	black	125052	4001363002041
	ocean	127741	4001363000177
	white (ET176)	125050	4001363002027
Dimensions	Ø 15/11 x 60 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	111330	4001363000191
	Skinfoam blue	111328	4001363000184
	Skinfoam violet (ET176)	128614	4001363003406
Material Colour	foam + imitation leather		
	jade	125055	4001363002065
	black	125056	4001363002072
	ocean	125053	4001363000207
	white (ET176)	125054	4001363002058
Dimensions	Ø 22/16,5 x 40 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115380	4001363000481
	Skinfoam blue	115381	4001363000498
	Skinfoam violet (ET176)	128615	4001363003413
Material Colour	foam + imitation leather		
	jade	125058	4001363002102
	black	103019	4001363000214
	ocean	127745	4001363002089
	white (ET176)	125057	4001363002096
Dimensions	Ø 22/16,5 x 60 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115382	4001363000504
	Skinfoam blue	115171	4001363001044
	Skinfoam violet (ET176)	128616	4001363003420
Material Colour	foam + imitation leather		
	jade	125061	4001363002126
	black	125062	4001363002133
	ocean	125059	4001363000221
	white (ET176)	125060	4001363002119



Trade purchase prices net per piece plus VAT


POSITIONING AIDS

Head supports



Occipital support, straight for adults

Dimensions	25 x 24,5 x 9,5 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	103869	4001363000245
	Skinfoam blue	103870	4001363000238
	Skinfoam violet (ET176)	108005	4001363003017
Material Colour	PE-thermoforming blue	100092	4001363000252



Occipital support, edges rounded for adults

Dimensions	23 x 23,5 x 9 cm	Item No.:	UDI
Material Colour	Soft-PE violet	114997	4001363001556

edges rounded



Occipital support, edges rounded for adults

Dimensions	21,5 x 23 x 9,5 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	107232	4001363001051
	Skinfoam blue	103874	4001363001068
	Skinfoam violet (ET176)	128617	4001363003437



Occipital support Polyprotect® for adults

Dimensions	24 x 25 x 10 cm	Item No.:	UDI
Material Colour	Polyprotect blue	103873	4001363000276



Occipital support Polyprotect® for adults for the operating theatre (viscoelastic/conductive)


Dimensions	24 x 25 x 10 cm	Item No.:	UDI
Material Colour	Polyprotect black	103872	4001363002652



Trade purchase prices net per piece plus VAT

Headrest for adults
left and right positioning with access for tubes

Dimensions	40 x 32 x 10 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115745	4001363000528
	Skinfoam blue	103876	4001363001082
	Skinfoam violet (ET176)	128618	4001363003444




Head support / abdominal Polyprotect®
for adults

Dimensions	22,5 x 22 x 11,5 cm	Item No.:	UDI
Material Colour	Polyprotect blue	115068	4001363002669



Occipital support, straight
for children

Dimensions	21,5 x 23 x 7 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	114566	4001363002645
	Skinfoam blue	107583	4001363000269
	Skinfoam violet (ET176)	128619	4001363003451



Occipital support, edges rounded
for children

Dimensions	21,5 x 23 x 7 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115744	4001363000511
	Skinfoam blue	103875	4001363001075
	Skinfoam violet (ET176)	128621	4001363003468



Trade purchase prices net per piece plus VAT

POSITIONING AIDS

Supporting aids & mattresses




Supporting aid supine position

Dimensions	40 x 25 x 12 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	108350 109630 108072	4001363001099 4001363001105 4001363003161
Dimensions	53 x 34 x 14 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	103821 115751 108074	4001363221112 4001363000740 4001363003178



Supporting aid lateral position

Dimensions	60 x 30 x 13 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	122487 109632 128622	4001363002676 4001363001129 4001363003475



Mattress

Dimensions	180 x 60 x 3 cm	Item No.:	UDI
Material Colour	imitation leather jade black	125090 125088	4001363002478 4001363002454
Plastic zip	ocean white (ET176)	125091 125089	4001363002485 4001363002461
Dimensions	200 x 60 x 3 cm	Item No.:	UDI
Material Colour	imitation leather jade black	125085 113945	4001363002508 4001363002492
Plastic zip	ocean white (ET176)	125087 125086	4001363002522 4001363002515



Trade purchase prices net per piece plus VAT

Arm pad with indentation

Dimensions	30 x 15 x 12/4 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	113760	4001363001136
	Skinfoam blue	115773	4001363000573
	Skinfoam violet (ET176)	128623	4001363003482
Material Colour	foam + imitation leather		
	jade	127653	4001363001891
	black	125063	4001363001907
	ocean	127647	4001363001808
	white (ET176)	127652	4001363001884



Arm pad with indentation

Dimensions	45 x 15 x 13/5 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	113762	4001363001143
	Skinfoam blue	103827	4001363001150
	Skinfoam violet (ET176)	128624	4001363003499
Material Colour	foam + imitation leather		
	jade	127656	4001363001839
	black	127659	4001363001877
	ocean	127655	4001363001853
	white (ET176)	127658	4001363001860



Arm pad with indentation

Dimensions	60 x 15 x 13/5 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115775	4001363000566
	Skinfoam blue	115776	4001363000580
	Skinfoam violet (ET176)	128625	4001363003505
Material Colour	foam + imitation leather		
	jade	125066	4001363001938
	black	125067	4001363001945
	ocean	127733	4001363001914
	white (ET176)	125065	4001363001921



Arm pad without indentation

Dimensions	30 x 15 x 12/4 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103829	4001363000009
Material Colour	Skinfoam yellow	127948	4001363002683
	Skinfoam blue	110053	4001363002690
	Skinfoam violet (ET176)	128627	4001363003512
Material Colour	foam + imitation leather		
	jade	125069	4001363001976
	black	125070	4001363001983
	ocean	125068	4001363001952
	white (ET176)	127735	4001363001969



POSITIONING AIDS

Triangles



Triangle

Dimensions	60 x 25 x 12,5 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	112193	4001363001167
	Skinfoam blue	107226	4001363001174
	Skinfoam violet (ET176)	108062	4001363003123



Triangle

Dimensions	36 x 24 x 12 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	112191	4001363001181
	Skinfoam blue	107224	4001363001198
	Skinfoam violet (ET176)	108064	4001363003130



Triangle

Dimensions	24 x 18 x 12 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103833	4001363001204
Material Colour	Skinfoam yellow	103835	4001363001211
	Skinfoam blue	103836	4001363001228
	Skinfoam violet (ET176)	108067	4001363003147
Material Colour	Polyprotect® blue	110038	4001363001822



Triangle

Dimensions	18 x 10 x 5 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103838	4001363001235
Material Colour	Skinfoam yellow	103840	4001363001242
	Skinfoam blue	103841	4001363001259
	Skinfoam violet (ET176)	108068	4001363003154



Triangle

Dimensions	24 x 18 x 12 cm	Item No.:	UDI
Material Colour	Soft-PE violet	122094	4001363001525


edges rounded



Trade purchase prices net per piece plus VAT

Rectangle

Dimensions	36 x 24 x 3 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115783	4001363000559
	Skinfoam blue	110654	4001363001266
	Skinfoam violet (ET176)	108018	4001363003024




Rectangle

Dimensions	36 x 24 x 5 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103845	4001363001273
Material Colour	Skinfoam yellow	110976	4001363001280
	Skinfoam blue	110036	4001363001297
	Skinfoam violet (ET176)	103846	4001363003000



Rectangle

Dimensions	18 x 24 x 3 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115786	4001363000597
	Skinfoam blue	110008	4001363001303
	Skinfoam violet (ET176)	108031	4001363003048



Rectangle

Dimensions	18 x 24 x 5 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103849	4001363001310
Material Colour	Skinfoam yellow	103850	4001363001327
	Skinfoam blue	103851	4001363001334
	Skinfoam violet (ET176)	108053	4001363003086
Material Colour	Polyprotect® blue	115201	4001363001785



Rectangle

Dimensions	18 x 24 x 5 cm	Item No.:	UDI
Material Colour	Soft-PE violet	122093	4001363001518

edges rounded



Trade purchase prices net per piece plus VAT

POSITIONING AIDS

Wedges



Wedge

Dimensions	15° 36 x 24 x 7/1 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	107809 110940 108028	4001363001341 4001363001358 4001363003031




Wedge

Dimensions	15° 18 x 28 x 8/1 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103855	4001363001365
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115795 103856 108057	4001363000603 4001363001372 4001363003093



Wedge

Dimensions	15° 18 x 24 x 7/1 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103858	4001363001389
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	110979 103859 108033	4001363001400 4001363001396 4001363003055
Material Colour	Polyprotect® hellblue	115797	4001363000610



Wedge


Dimensions	20° 18 x 28 x 10/1 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	115799 110651 108044	4001363000627 4001363001419 4001363003062



Wedge

Dimensions	15° 18 x 24 x 7/1 cm	Item No.:	UDI
Material Colour	Soft-PE violet	114999	4001363001532


edges rounded



Trade purchase prices net per piece plus VAT

Circular foam pad

Dimensions	Ø15 x 4 cm	Item No.:	UDI
Material Colour	Skinfoam yellow	115807	4001363000535
	Skinfoam blue	115815	4001363000542
	Skinfoam violet (ET176)	108058	4001363003109




Circular foam pad

Dimensions	Ø22 x 5 cm	Item No.:	UDI
Material Colour	PE-thermoforming blue	103861	4001363001426
	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	103862	4001363001433
103863		4001363001440	
108060		4001363003116	



Circular foam pad with hole

Dimensions	Ø22 x 5 cm	Item No.:	UDI
Material Colour	Skinfoam yellow Skinfoam blue Skinfoam violet (ET176)	103864	4001363001457
		103865	4001363001464
		108051	4001363003079



Trade purchase prices net per piece plus VAT

POSITIONING AIDS

Positioning-aids-set



Positioning-aids-set | 10 pieces in a plastic case

Material Colour	PE-thermoforming blue	109633	4001363002706
Material Colour	Skinfoam yellow Skinfoam blue (ET176)	110913 113577	4001363000030 4001363000047

Item 113577 – Skinfoam blue
consists of following items:

1 x 103863	circular foam pad	22 x 5 cm
1 x 103865	circular foam pad with hole	22 x 5 cm
1 x 103893	roll	10 x 25 cm
1 x 103851	rectangle	18 x 24 x 5 cm
2 x 103836	triangle 45°	18 x 24 x 12 cm
2 x 103841	triangle 45°	18 x 10 x 5 cm
2 x 103859	wedge 15°	18 x 24 x 7/1 cm
in a plastic case		

Item 109633 – PE-thermoforming, blue
consists of following items:

1 x 103829	arm pad	30 x 15 x 12/4 cm
1 x 103845	rectangle	36 x 24 x 5 cm
1 x 103855	wedge 15°	18 x 28 x 8/1 cm
1 x 103858	wedge 15°	24 x 18 x 7/1 cm
2 x 103833	triangle 45°	24 x 18 x 12 cm
2 x 103838	triangle 45°	18 x 10 x 5 cm
2 x 103849	rectangle	24 x 18 x 5 cm
in a plastic case		



Trade purchase prices net per piece plus VAT

Positioning-aids-set | 11 pieces
in a plastic case

Material Colour	PE-thermoforming blue	109034	4001363001471
Material Colour	Skinfoam yellow Skinfoam blue (ET176)	114841 109631	4001363001488 4001363001495

Item 109631 - Skinfoam blue consists of following items:

1 x 103863	circular foam pad	22 x 5 cm
1 x 103865	circular foam pad with hole	22 x 5 cm
1 x 103893	roll	10 x 25 cm
1 x 103851	rectangle	18 x 24 x 5 cm
2 x 103836	triangle 45°	18 x 24 x 12 cm
2 x 103841	triangle 45°	18 x 10 x 5 cm
2 x 103859	wedge 15°	18 x 24 x 7/1 cm
1 x 103870	occipital support	25 x 24,5 x 9,5 cm

in a plastic case

Item 109034 - PE-thermoforming, blue consists of following items:

1 x 103829	arm pad	30 x 15 x 12/4 cm
1 x 103845	rectangle	36 x 24 x 5 cm
1 x 103855	wedge 15°	18 x 28 x 8/1 cm
1 x 103858	wedge 15°	24 x 18 x 7/1 cm
2 x 103833	triangle 45°	24 x 18 x 12 cm
2 x 103838	triangle 45°	18 x 10 x 5 cm
2 x 103849	rectangle	24 x 18 x 5 cm
1 x 100092	occipital support	25 x 24,5 x 9,5 cm

in a plastic case



Trade purchase prices net per piece plus VAT

POSITIONING AIDS

Positioning-aids-set



Positioning-aids-set | 8 pieces in a plastic case

Material Colour	Soft-PE violet	122770	4001363001501
--------------------	-------------------	--------	---------------

Item 122770 - Soft-PE violet
consists of following items:

1 x 122092	1/2 roll	Ø 15/7,5 x 40 cm
2 x 122093	rectangle	18 x 24 x 5 cm
2 x 122094	triangle 45°	24 x 18 x 12 cm
2 x 114999	wedge 15°	18 x 24 x 7/1 cm
1 x 114997	occipital support - rounded	25 x 24,5 x 9,5 cm

in a plastic case



Trade purchase prices net per piece plus VAT


Scrotal support | 4 pieces

Dimensions	13 x 11 x 8 cm	Item No.:	UDI
Material Colour	ET127 blue	132248	4001363004939




Foot and elbow protection

		Item No.:	UDI
Material Colour	ET127 blue	132249	4001363004946




Wedge for 30°- positioning

Dimensions	60 x 30/12 x 12 cm	Item No.:	UDI
Material Colour	ET127 blue	103811	4001363004984



Topper - profiled foam

Material Colour	ET127 blue	Item No.:	UDI
Dimensions	60 x 60 x 5 cm	103922	4001363004953
Dimensions	90 x 60 x 5 cm	103923	4001363004968
Dimensions	195 x 90 x 5 cm	103924	4001363004977



Heel support | pair

Material Colour	ET127 blue	Item No.:	UDI
Dimensions	10/3 x 4 cm	132250	4001363004991
Dimensions	15/5 x 4 cm	132251	4001363005004
Dimensions	20/8 x 4 cm	132252	4001363005011



Trade purchase prices net per piece plus VAT

TRAINING ARTICLES

applicable for sports,
leisure and training activities



Polyform-Relax roll

Dimensions	Ø 6 x 9 cm	Item No.:	
Material Colour	Relax 36 very soft, yellow	103883	
Material Colour	Relax 38 soft, pink	103884	
Material Colour	Relax 41 middle, blue	103885	Price: 4,01 EUR



Handballs | 10 pieces

Dimensions	Ø 7 cm	Item No.:	
Material Colour	Polyurethane yellow	103886	




Handrolls | 10 pieces

Dimensions	Ø 6 x 9 cm	Item No.:	
Material Colour	Polyurethane yellow	103887	



Handtrainer | 10 pieces

Material Colour	Polyethylene mixed (red+blue)	112080	
Material Colour	Polyethylene middle (red)	103888	
Material Colour	Polyethylene hard (blue)	103889	



Handle thickenings Polyethylen

Quantity Colour	10 pieces white	103890	
Quantity Colour	25 pieces white	115824	



Trade purchase prices net per piece plus VAT

No medical products according to Regulation (EU) 2017/745

Cotton, white



Imitation leather, white



Imitation leather, black




Imitation leather, jade



Imitation leather, ocean



Skinfoam, yellow



Skinfoam, blue




Skinfoam, violet




CUSTOM-MADE PRODUCTS

Individual form pads (from 1 piece)


Type: F-Wedge




Type: Z-Wedge




Type: 1/2 roll




Type: H-Triangle



Type: U-Triangle



Type: Rectangle



Please specify the desired type, the corresponding dimensions, the foam quality, the surface and the quantity.

COVERING FABRICS & COATINGS

Cotton / Imitation leather / Skinfoam



Description

Skinfoam is a plastic coating which is applied thinly to the foam, to create a disinfectable surface, which can be wiped with a damp cloth. The surface is not waterproof and has, depending on stress, a limited lifespan.

Technically, all foam qualities can be coated.

We do not recommend the coating of composite foams, though, as the coating looks poorly on foam qualities with an „uneven“ surface. Furthermore, the coating of viscoelastic foams is not recommended, due to the fact that they „break“ quite quickly at the edges / corners.


Ordinary polyurethane or polyethylene foams can be coated without any problems. However, it should be mentioned that the coating of polyethylene surfaces is more expensive because the surfaces have to be pretreated.

Prices by expenses! Please inquire separately the surfaces to be coated, stating size (sketch) and the corresponded foam quality (PE or PU).

Standard colours: yellow, blue (radiolucent), violet

Care-recommendation:

Skinfoam-coating
(Made of Polyurethane)



advantages:

- > skin-friendly Polyurethane-coating
- > durable surface
- > custom-made products are possible

X-ray possibility:

- > from 70 –140 KV
- > colours: yellow & blue

colours:

- > usable for x-ray

care:

- > wipeable with a damp cloth


disinfection:

undiluted Antifect Liquid (Schülke & Mayr)

ATTENTION:

The material is moisture-repellent but not absolutely waterproof. Please avoid waterlogging to prolong the life of your Polyform-product. We cannot guarantee a cleaning of some substances(e.g. blood, secretion etc.).





POLYFORM Kunststofftechnik
GmbH & Co. Betriebs KG
Braasstraße 15, D- 31737 Rinteln
Tel.: +49 (0) 57 51 / 40 08-600
Fax: +49 (0) 57 51 / 40 08-699
med.formpolster@polyform.de
www.polyform.de